



## CELOX<sup>™</sup> Z-Fold Gauze (5ft)

### How to Use/Recommendations for Removal




- 1** Tear open pack, take one end of the CELOX<sup>™</sup> Z-Fold Gauze while applying pressure to the wound.



- 2** Tightly pack the unfolding CELOX<sup>™</sup> Z-FOLD Gauze directly to the source of bleeding, feeding in the CELOX<sup>™</sup> gauze until the wound is filled above skin level. Excess gauze can be cut or torn.



- 3** Apply FIRM pressure for 3 minutes or until bleeding stops.



- 4** Secure with a dressing to maintain pressure and discard any remaining CELOX<sup>™</sup> Z-Fold Gauze.

- 5** If not in a medical facility, transfer patient as soon as possible and show empty packet to medical personnel.

#### ATTENTION MEDICAL PERSONNEL:

- 1.** Physically remove CELOX<sup>™</sup> Z-Fold Gauze from the wound and any loose surface granules.
- 2.** Fully flood entire wound area with sterile saline irrigation solution.
- 3.** Proceed with normal irrigation and/or suction.
- 4.** Ensure all product is removed from the wound prior to initiation of wound treatment.

**Warnings:** Device is not intended for surgical use. Do not use in abdominal wounds and wounds not amenable to pressure. Loss of sterility potentially poses a risk of infection. Do not apply over eyes. If eye irritation occurs flush with water for 5 minutes. Contains chitosan from shellfish. Allergy studies show no adverse reaction to date. Data on file at Medtrade Products Ltd.

## PRODUCT EVALUATION FORM

Evaluation forms can be emailed to [info@omni-stat.com](mailto:info@omni-stat.com) or given to your local sales representative

Facility/Hospital Name: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_ Provider Signature: \_\_\_\_\_

### Product Preparation Used

- OMNI-STAT Granules (3g)
- CELOX™ Z Fold Gauze (5ft)
- CELOX™ Granules (15g)
- OMNI-STAT Gauze (4in x 4in)
- CELOX™-A (6g)
- CELOX™ RAPID


### Type of Wound

- Gunshot
- Laceration/Abrasion
- Puncture
- Blunt Trauma
- Skin Tear
- Diabetic Foot Ulcer
- Pressure Ulcer
- Other: \_\_\_\_\_

### Hold Pressure Time

- 1 minute
- 3 minutes
- 5 minutes
- 10 minutes
- Other: \_\_\_\_\_

### Wound Location



Please circle the location of the wound on the diagram above

### Was the Patient Anticoagulated?

- YES,  
If Yes, please list type  
of anticoagulant: \_\_\_\_\_
- NO

---

---

---

### Was hemostasis achieved?

- YES
- NO

---